

Specifications

Visit our website at easystand.com for complete specifications on all EasyStand products and options.



	Bantam Extra Small	Bantam Small	Bantam Medium
Approximate Height Range	28"-40" (71-102cm)	36"-54" (91-137cm)	4'-5'6" (122-168cm)
Maximum User Weight	50 lbs (23kg)	100 lbs (45kg)	200 lbs (91kg)
Seat to Foot Plate Range	4.5"-16.5" (11-42cm)	4.5"-16.5" (11-42cm)	11"-19.5" (28-50cm)
Seat Depth Range w/Back	7"-12" (18-31cm)	11"-16" (28-41cm)	14"-19" (36-48cm)
Seat Height from Floor	19" (48cm)	19" (48cm)	21.5" (55cm)
Knee Pad Size	2.5"Wx5"H (6x13cm)	3.5"Wx5"H (9x13cm)	7"Wx9"H (18x23cm)
Knee Pad Depth Range	2.5"-8" (6-20cm)	2.5"-8" (6-20cm)	3"-6" (8-15cm)
Weight of Basic Unit	54 lbs (24.5kg)	54 lbs (24.5kg)	108 lbs (49kg)
Frame Footprint	24.5"x36.5" (62x93cm)	24.5"x36.5" (62x93cm)	24.5"x44.5" (61x113cm)
Supine	0-90°	0-90°	0-90°

	Evolv Medium	Evolv Large	Evolv XT	Glider Medium	Glider Large	StrapStand
Approximate Height Range	4'-5'6" (122-168cm)	5'-6'2" (152-188cm)	6'-6'10" (183-211cm)	4'-5'6" (122-168cm)	5'-6'2" (152-188cm)	5'-6'5" (152-195cm)
Maximum User Weight	200 lbs (90kg)	280 lbs (127kg)	350 lbs (159kg)	200 lbs (90kg)	280 lbs (127kg)	350 lbs (159kg)
Seat to Foot Plate Range	11.5"-18" (29-46cm)	14"-20.5" (36-52cm)	17"-23" (43-58cm)	13.5"-16.5" (34-42cm)	15.5"-18.5" (39-47cm)	N/A
Seat Depth Range w/Back	14"-19" (36-48cm)	18"-23" (46-59cm)	20"-25" (50-64cm)	14"-19" (36-48cm)	18"-23" (46-58cm)	N/A
Seat Height from Floor	21.5" (55cm)	21.5" (55cm)	23.5" (60cm)	21.5" (55cm)	21.5" (55cm)	N/A
Knee Pad Size	17"Wx9"H (43x23cm)	17"Wx9"H (43x23cm)	17"Wx9"H (43x23cm)	7"Wx7"H (18x18cm)	7"Wx7"H (18x18cm)	17"Wx9"H (43x23cm)
Knee Pad Depth Range	2"-7" (5-18cm)	2"-7" (5-18cm)	2.5"-8" (6-20cm)	N/A	N/A	6"-9.5" (15-24cm)
Weight of Basic Unit	101 lbs (46kg)	104 lbs (47kg)	112 lbs (51kg)	174 lbs (79kg)	174 lbs (79kg)	125 lbs (58kg)
Frame Footprint	26.5"x36" (67x91cm)	26.5"x36" (67x91cm)	28.5"x38" (72x97cm)	26.5"x41" (67x104cm)	26.5"x41" (67x104cm)	28"x34" (71x86cm)